SOFTWARE



Open Access



CellProfiler Tracer: exploring and validating high-throughput, time-lapse microscopy image data

Mark-Anthony Bray and Anne E. Carpenter*

Abstract

Background: Time-lapse analysis of cellular images is an important and growing need in biology. Algorithms for cell tracking are widely available; what researchers have been missing is a single open-source software package to visualize standard tracking output (from software like CellProfiler) in a way that allows convenient assessment of track quality, especially for researchers tuning tracking parameters for high-content time-lapse experiments. This makes quality assessment and algorithm adjustment a substantial challenge, particularly when dealing with hundreds of time-lapse movies collected in a high-throughput manner.

Results: We present CellProfiler Tracer, a free and open-source tool that complements the object tracking functionality of the CellProfiler biological image analysis package. Tracer allows multi-parametric morphological data to be visualized on object tracks, providing visualizations that have already been validated within the scientific community for time-lapse experiments, and combining them with simple graph-based measures for highlighting possible tracking artifacts.

Conclusions: CellProfiler Tracer is a useful, free tool for inspection and quality control of object tracking data, available from http://www.cellprofiler.org/tracer/.

Keywords: CellProfiler, Time lapse, Quality assessment, Fluorescence microscopy, Image analysis, Data visualization, Data exploration

Background

Time-lapse assays probe biological questions that can only be investigated by observing the dynamic behavior of organisms, cells, organelles, or molecular assemblies over time [1]. The combination of automated imaging and large-scale, high-content, live-cell experiments is capable of delivering large amounts of data in very little time [2]. However, time-lapse imaging is acutely susceptible to many artifacts that negatively affect the proper identification and tracking of cells; the appearance of such anomalies in a single frame can ruin an entire time series. Thus, image and image analysis quality requirements for time-lapse microscopy are more stringent and, due to the volume of data, automated quality control is more necessary.

Interfaces for review or correction of time-lapse data are sometimes provided in customizable open-source software but are usually manual in nature, requiring

* Correspondence: anne@broadinstitute.org

visual inspection to detect aberrations [3, 4]. Alternately, commercial software may present such functionality within a polished interface (e.g., Imaris by Bitplane, Volocity by Perkin-Elmer, MetaMorph by Molecular Devices), but such packages are not open-source, precluding access to or adjustment of features and underlying algorithms. We saw a need for a tool that would link cell images themselves directly to their morphological measurements within a tracking assessment tool; without this, a valuable opportunity is missed for the researcher to visually assess important changes in cell morphology and cellular context that accompany particular tracking results.

The CellProfiler biological image analysis package is widely used for collecting an extensive suite of morphological, intensity, and textural features for cells and organisms in high-content screens [5, 6]. Moreover, it includes modular cell-tracking capabilities for time-lapse assays, such as the linear assignment problem (LAP) approach [7], which provides robust tracking by closing temporal gaps and capturing object merges and splits.



© 2015 Bray and Carpenter. **Open Access** This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

Broad Institute of MIT and Harvard, 415 Main St, Cambridge, MA 02142, USA

CellProfiler is one of the few options for conveniently combining the need for robust cellular identification and the ability to process large numbers of time-lapse movies [8-10]. However, configuring tracking algorithm parameters has been tedious without a tool to readily assess track quality. We present CellProfiler Tracer to enable the visualization of the rich set of cellular features characteristic of high-content time-lapse assays, as well as to provide measures for assessing track quality.

Implementation

CellProfiler Tracer is implemented as part of the CellProfiler Analyst software package [11] (available from http:// www.cellprofiler.org). Although most seamlessly used with data from CellProfiler, the Tracer software is intended as a visualization and quality assessment tool compatible with high-content object tracking data derived from any twodimensional time-lapse image sets, which most commonly involve fluorescence or brightfield microscopy. Thus, Tracer is not itself a tracking algorithm, nor a generalpurpose image visualization tool, nor a tool for manual track editing, but may be used in conjunction with other software for those purposes [12, 13]. We used CellProfiler Analyst as the foundation for Tracer, as CellProfiler Analyst was designed for visualizing large, multi-parametric data sets, with the ability to create various plots of cellular features. It emphasizes linking the plotted data to the originating image for visual inspection and improved biological interpretation. To create Tracer, we added visualizations to CellProfiler Analyst that are specific to time-lapse tracking data and that have already been proven useful in the scientific community but have not as yet existed in a single freely-available and open-source software package, namely the following:

- XYT plot (Fig. 1a): An *XYT plot* is a 3-D plot of the cell centroid versus time and is useful as a straight-forward means of visualizing discrete cellular trajectories [14, 15]. Tracer can color-code the trajectories based on a selected per-cell feature, so that variations in the trajectory color reflect the size, shape, intensity, or other high-content, multi-channel features collected during the experiment. The plot can be rotated using the mouse to view the trajectories from any angle.
- Lineage tree (Fig. 1b): *Lineage trees* omit the positional information of XYT plots to display relationships between cell descendants over time; such graphs are commonly used for developmental mechanism and cell cycle



here by the frame number; the trajectories can be color-coded to be any cell measurement of interest; (**b**) the lineage tree panel, highlighting the ancestor/progeny relationships corresponding to the trajectories in (**a**), and (**c**) the control panel containing various display tools. Other visualizations include (**d**) synchrograms of selected cells, as well as heatmaps (shown in Fig. 2) progression studies [4, 16–19]. Each cell at a given timepoint is represented as a node, connected by edges to the tracked predecessors and successors. As with the XYT plot, the nodes can be color-coded according to the desired percell image feature. Errors in object segmentation often appear and disappear within a few frames, which results in an object split followed by a remerge (or vice versa). This topological distortion is most obvious in the lineage display and is used to call attention to the need for adjustment to the image analysis parameters.

- Synchrograms (Fig. 1d): A *synchrogram* is a sequence of images of an individual cell over time. This visualization helps track the progression of subcellular processes [20]. The selected cell is centered in each frame to remove motion as a visual degree of freedom; this also makes frame-to-frame tracking errors immediately apparent. This tool allows for follow-up on suspect trajectories or lineages identified with the other tools, using a simple point-and-click interface linked to the other plots.
- Heatmaps (Fig. 2): *Heatmaps* are used to represent numerical data graphically as a colored two-dimensional matrix. Tracer can average each perobject feature across all trajectories at each timepoint and display a heatmap of the result. This display provides a simple means for the user to visually evaluate the data for significant trends in the cell population, even if the time-lapse data is not temporally synchronous. This can be helpful for quality control purposes (e.g., a given timepoint was transiently out of focus) or perhaps relevant to the phenotype under investigation (e.g., observing trends in response to drug treatment).

A control panel (Fig. 1c) is used to adjust various aspects of the display panels. Other visualization features, such as plots and trajectory selections, are provided via a context menu available on the display panels (Additional file 1: Figure S1). In addition, diagnostics related to the LAP tracking method [21] and other quality metrics (described below) may also be displayed as an aid to tracking optimization (Additional file 1: Figure S2).

Tracer is designed to process and explore any MySQLor SQLite-based database of image-based screening data structured according to the following simple schema:

• An image table where each row corresponds to an image acquired at a unique timepoint and field of view and the columns contain the image data (e.g., the name of the treatment condition, the path to and filename of the original image, etc.). A



Fig. 2 Heatmap of high-content cellular time-lapse measurements. The per-nucleus measurements from a *Drosophila* time-lapse movie are averaged over all nuclei for each timepoint; the measurements were collected by CellProfiler software. Feature values were normalized from 0 to 1 for visualization purposes. Feature names were omitted for conciseness but are provided in the Tracer display; the features shown are listed in order in the Additional file 3: Table S1, and are further described in the CellProfiler documentation

requirement for this table is an image index, given as a column of integers referencing each site (i.e., field of view) acquired.

• An object table in which each row represents an object (e.g., cells) from a given image and the columns contain the collected object measurements (e.g., area of the cell, intensity of DNA stain in the nucleus, location of the cell in the original image). Required for this table is an image index as described above, as well as an object index given as a column of integers referencing each object identified in an image. An (x,y) location for each object is also required (e.g., the cell centroid) as columns in the table; this permits limited 3D + *t* functionality if a 2D centroid of each 3D object is

provided rather than an entire 3D segmentation, e.g., by using a maximum projection into the XY plane or a particular Z-slice.

• In addition, an object relationship table is also required. Each row corresponds to the image and object index for a given object and that of its "parent", i.e., tracked predecessor.

The above data tables can be automatically generated by CellProfiler using its ExportToDatabase and ExportToSpreadsheet modules; example files are provided at http:// www.cellprofiler.org/tracer/. However, these schema are intended to be simple enough that a third-party software package (e.g., MATLAB, another object tracking package) can easily format its data accordingly. Using this format, a trajectory (defined as the frame-to-frame path followed by an object over time) can be captured by following the ancestor-progeny mapping for an object. The relationship table captures both one-to-many and many-to-one object mappings created by splits and merges, as well as temporal gaps produced by transient object disappearances.

The primary challenge in handling aberrations in time-lapse data is discriminating between natural biological behaviors versus analysis artifacts (e.g., cell division during mitosis versus improper cell splitting due to mis-segmentation). To address this issue, Tracer allows the user to assess the quality of object trajectories by treating the trajectories as a network graph and highlighting possible aberrations in the graph connectivity. The use of graph-based methods to analyze and resolve defects in time-lapse tracking data has been explored previously [16, 22–25]; here, we use Tracer to simply bring attention to aspects of the network that may represent possible mistakes in cellular segmentation and tracking.

Results and discussion

We demonstrate the data analysis and quality assessment features of CellProfiler Tracer using two time-lapse movies that reveal the dynamics of nuclear division: a Drosophila blastoderm embryo with GFP-histone marking the nuclear DNA (Foe lab, University of Washington, unpublished data used with permission) and MCF-7 nuclei tagged with NLS-mCerulean fusion protein (Ramaswamy lab, Massachusetts General Hospital Cancer Center, unpublished data used with permission); further details on these data and permissions for use are included in the Additional file 2. We used CellProfiler to identify the nuclei, track them over time, and measure over 130 features of area, shape, intensity and texture (see the Additional file 3: Table S1 for the full list of cellular features); the image data and CellProfiler pipelines are available from http:// www.cellprofiler.org/tracer/. In the case of the Drosophila embryos, the nuclei proceed through the cell cycle in synchrony, due to their sharing a common cytoplasm. From the heatmap shown in Fig. 2, one can see that the rich set of features derived from GFP-histone expression at each time point of the movie could be used to fingerprint nuclei at particular phases of the cell cycle. Similarly, any of the features collected can be visualized on the XYT and lineage panels by selecting the desired measurement for color-coding the object tracks; Additional file 1: Figure S3 shows examples for the two data sets. Hypotheses about cell behavior can be generated and tested using this view.

Cells typically exhibit a limited range of dynamic behavior, and hence the resultant network graphs are expected to assume only certain topologies, as shown in a movie of MCF-7 nuclei (Fig. 3a). Therefore, deviations from expected topologies (Fig. 3b) may indicate that something is amiss in the tracking. The user can select to display three different graph deviations in the Tracer interface: loops, crossings, and singletons; all of these are evaluated when the data is first loaded. For example, transient split/merge ("loops") or merge/split ("crossings") events are unlikely to occur in typical biological settings and may indicate an object mis-segmentation (Fig. 3c, d). Likewise, a very short trajectory ("singleton", although the precise number of frames can be selected by the researcher) may correspond to a spurious object detection. For singletons, the user can produce a display showing the distribution of track lengths and a chart listing the total number of tracks and the median, 10th and 90th percentiles of the track lengths (Additional file 1: Figure S4).

CellProfiler Tracer examines the tracking data for these deviations and highlights suspect nodes by colorcoding them in the display panels. In the case of singletons, these nodes may be removed to create a new network graph saved with the original data for retrieval across Tracer sessions. It should be emphasized that the highlighted nodes are intended to call attention to possible object segmentation or tracking problems; the absence of suspect nodes in a particular dataset does not preclude other errors. However, a large number of highlighted nodes would indicate the need for further optimization of the original cell segmentation settings. While other tracking interfaces include manual or semisupervised segmentation editing [4, 26], this functionality is currently outside the scope of Tracer. This is because the tool aims to support the completely automated analysis of thousands of time-lapse image sets, for which manual intervention is impossible and some amount of automated error must be tolerated. Thus, the tool is designed to assist researchers in selecting appropriate segmentation and tracking parameters in other software (like CellProfiler), such that the resulting data files are as high quality as possible. If feasible and necessary (e.g., for small-scale experiments), manual editing of individual trajectories can be carried out using other software, after



Tracer has assisted in optimizing the automated segmentation and tracking parameters.

To illustrate the use of CellProfiler Tracer for optimizing parameters, we tracked nuclei in the MCF-7 timelapse data set using CellProfiler; for this movie, the default settings for its LAP tracking method yielded substantial errors (Additional file 1: Figure S5). Anecdotally, researchers have reported that the lack of visual feedback makes adjusting the algorithm's many parameters extremely challenging. Using Tracer to visualize and assess the resulting tracks, the tracking settings in Cell-Profiler were then changed accordingly and the nuclei re-analyzed to produce revised tracks. This procedure was iteratively repeated until the final tracks reflected the actual temporal behavior of the nuclei; these optimized settings were then confirmed in Tracer to reproduce the expected tracking behavior for a different MCF-7 data set (Additional file 1: Figure S6). It bears noting that while only the lineage panel is shown in Additional file 1: Figure S5 for brevity, all three visualization tools were employed for optimization. While this one particular movie could have been analyzed more quickly using a tool for manual correction; the value in using Tracer is to optimize automated settings on subsets of time-lapse data, so that they can be applied to hundreds of thousands of data sets, for example, using CellProfiler's high-throughput interface.

Conclusions

As cellular tracking matures and the size of microscopy data sets continues to increase, progress in validating tracking quality will make powerful time-lapse experiments on larger data sets feasible [26–29]. CellProfiler Tracer is a tool that augments the cell-tracking functionality of the CellProfiler biological image analysis package by visualizing multi-parametric time-lapse data. The software incorporates graph-based assessment of tracking quality, and makes it easy to produce and interact with XYT plots, lineage trees, synchrograms and heatmaps — visualizations that have proven useful but have not as yet existed in a single freely-available and opensource software package.

Availability and requirements

- Project name: CellProfiler Tracer
- Project home page: http://www.cellprofiler.org/ tracer/ (installer), https://github.com/CellProfiler/ CellProfiler-Analyst/tree/cellprofiler-tracer (source code)
- Operating systems: Windows, 64-bit

- Other requirements: The Tracer source code for CellProfiler Analyst requires the following libraries (see the project page for the most up-to-date requirements):
 - For basic CellProfiler Analyst functionality
 - Python 2.8 or greater (3.0 is not currently supported)
 - NumPy 1.71 or greater
 - SciPy
 - wxPython
 - scikit-learn
 - MySOLdb
 - matplotlib
 - javabridge
 - python-bioformats
 - verlib (required by distutils)
 - \bigcirc For Tracer functionality
 - Enthought Tool Suite (for Mayavi2)
 - VTK, 5.10 or greater
 - NetworkX, 1.7 or greater
 - configobj (required by Enthought)
- License: GNU General Public License, Version 2.
- Any restrictions to use by non-academics: None

Additional files

Additional file 1: Figures S1–S6. Referenced in the main manuscript text. (PDF 742 kb)

Additional file 2: Details of the cellular image data including cell **type, acquisition technique, resolution and temporal resolution.** (PDF 119 kb)

Additional file 3: Listing of high-content image-based features generated by CellProfiler, categorized by feature type. (XLS 33 kb)

Competing interests

The authors declare that they have no competing interests.

Authors' contributions

MAB developed and implemented the software, and drafted the manuscript. AEC contributed to the design of the software and case studies and was involved in writing and editing the manuscript. Both authors reviewed and approved the final manuscript.

Acknowledgments

The authors thank the members of the Foe and Ramaswamy laboratories for contributing image data to the development of the software and this article, as well as John Albeck, Albert Yeh, and Evgeny Shlevkov for their helpful feedback on the Tracer tool. We are also grateful to Lee Kamentsky and Vebjørn Ljoså for helpful advice, coding edits and discussions on compilation, and the Imaging Platform team members for comments during manuscript preparation.

This work was supported by the National Science Foundation (NSF RIG DBI 1119830 to MAB) and National Institutes of Health (NIH R01 GM089652 to AEC).

Received: 8 April 2015 Accepted: 3 October 2015 Published online: 04 November 2015

References

 Spiller DG, Wood CD, Rand DA, White MRH. Measurement of single-cell dynamics. Nature. 2010;465:736–45.

- Neumann B, Walter T, Hériché J-K, Bulkescher J, Erfle H, Conrad C, et al. Phenotypic profiling of the human genome by time-lapse microscopy reveals cell division genes. Nature. 2010;464:721–7.
- Meijering E, Dzyubachyk O, Smal I. Methods for cell and particle tracking. Methods Enzymol. 2012;504:183–200.
- Winter M, Wait E, Roysam B, Goderie SK, Ali RAN, Kokovay E, et al. Vertebrate neural stem cell segmentation, tracking and lineaging with validation and editing. Nat Protoc. 2011;6:1942–52.
- Carpenter AE, Jones TR, Lamprecht MR, Clarke C, Kang IH, Friman O, et al. Cell Profiler: image analysis software for identifying and quantifying cell phenotypes. Genome Biol. 2006;7:R100.
- Kamentsky L, Jones TR, Fraser A, Bray M-A, Logan DJ, Madden KL, et al. Improved structure, function and compatibility for Cell Profiler: modular high-throughput image analysis software. Bioinformatics. 2011;27:1179–80.
- Jaqaman K, Loerke D, Mettlen M, Kuwata H, Grinstein S, Schmid SL, et al. Robust single-particle tracking in live-cell time-lapse sequences. Nat Methods. 2008;5:695–702.
- Failmezger H, Fröhlich H, Tresch A. Unsupervised automated high throughput phenotyping of RNAi time-lapse movies. BMC Bioinformatics. 2013;14:292.
- Rees P, Wills JW, Brown MR, Tonkin J, Holton MD, Hondow N, et al. Nanoparticle vesicle encoding for imaging and tracking cell populations. Nat Methods. 2014;11:1177–81.
- Frechin M, Stoeger T, Daetwyler S, Gehin C, Battich N, Damm E-M, et al. Cell-intrinsic adaptation of lipid composition to local crowding drives social behaviour. Nature. 2015;523:88–91.
- Jones TR, Kang IH, Wheeler DB, Lindquist RA, Papallo A, Sabatini DM, et al. Cell Profiler Analyst: data exploration and analysis software for complex image-based screens. BMC Bioinformatics. 2008;9:482.
- Chenouard N, Smal I, de Chaumont F, Maška M, Sbalzarini IF, Gong Y, et al. Objective comparison of particle tracking methods. Nat Methods. 2014;11:281–9.
- Peng H, Bria A, Zhou Z, lannello G, Long F. Extensible visualization and analysis for multidimensional images using Vaa3D. Nat Protoc. 2014;9:193–208.
- Li K, Miller ED, Chen M, Kanade T, Weiss LE, Campbell PG. Cell population tracking and lineage construction with spatiotemporal context. Med Image Anal. 2008;12:546–66.
- Tvaruskó W, Bentele M, Misteli T, Rudolf R, Kaether C, Spector DL, et al. Time-resolved analysis and visualization of dynamic processes in living cells. Proc Natl Acad Sci U S A. 1999;96:7950–5.
- Amat F, Lemon W, Mossing DP, McDole K, Wan Y, Branson K, et al. Fast, accurate reconstruction of cell lineages from large-scale fluorescence microscopy data. Nat Methods. 2014;11:951–8.
- Al-Kofahi O, Radke RJ, Goderie SK, Shen Q, Temple S, Roysam B. Automated cell lineage construction: a rapid method to analyze clonal development established with murine neural progenitor cells. Cell Cycle. 2006;5:327–35.
- Murray JI, Bao Z, Boyle TJ, Waterston RH. The lineaging of fluorescentlylabeled Caenorhabditis elegans embryos with StarryNite and AceTree. Nat Protoc. 2006;1:1468–76.
- Megason SG. In toto imaging of embryogenesis with confocal time-lapse microscopy. Methods Mol Biol. 2009;546:317–32.
- Sigal A, Milo R, Cohen A, Geva-Zatorsky N, Klein Y, Alaluf I, et al. Dynamic proteomics in individual human cells uncovers widespread cell-cycle dependence of nuclear proteins. Nat Methods. 2006;3:525–31.
- Jaqaman K, Danuser G. Computational image analysis of cellular dynamics: a case study based on particle tracking. Cold Spring Harb Protoc. 2009;2009:db.top65.
- Liu M, Yadav RK, Roy-Chowdhury A, Reddy GV. Automated tracking of stem cell lineages of Arabidopsis shoot apex using local graph matching. Plant J. 2010;62:135–47.
- Santella A, Du Z, Bao Z. A semi-local neighborhood-based framework for probabilistic cell lineage tracing. BMC Bioinformatics. 2014;15:217.
- Schiegg M, Hanslovsky P, Haubold C, Koethe U, Hufnagel L, Hamprecht FA. Graphical model for joint segmentation and tracking of multiple dividing cells. Bioinformatics. 2015;31:948–56.
- Schiegg M, Heuer B, Haubold C, Wolf S, Koethe U, Hamprecht FA. Proofreading guidance in cell tracking by sampling from tracking-by-assignment models. In 2015 IEEE 12th International Symposium on Biomedical Imaging (ISBI). IEEE:394–398.

- Wait E, Winter M, Bjornsson C, Kokovay E, Wang Y, Goderie S, et al. Visualization and correction of automated segmentation, tracking and lineaging from 5-D stem cell image sequences. BMC Bioinformatics. 2014;15:328.
- Maška M, Ulman V, Svoboda D, Matula P, Matula P, Ederra C, et al. A benchmark for comparison of cell tracking algorithms. Bioinformatics. 2014;30:1609–17.
- Chen Y, Ladi E, Herzmark P, Robey E, Roysam B. Automated 5-D analysis of cell migration and interaction in the thymic cortex from time-lapse sequences of 3-D multi-channel multi-photon images. J Immunol Methods. 2009;340:65–80.
- 29. Adanja I, Megalizzi V, Debeir O, Decaestecker C. A new method to address unmet needs for extracting individual cell migration features from a large number of cells embedded in 3D volumes. PLoS One. 2011;6:e22263.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar

BioMed Central

(

• Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit